

# Product Sheet

## H\_PDL1 MC38(mouse PDL1 KO) Cell Line

Catalog number: GM-C20064

Version 3.3.1.260512

<b>Description</b>	H_PDL1 MC38(mouse PDL1 KO) Cell Line is a clonal stable MC38 cell line that constitutively expresses human PDL1 gene, constructed using lentiviral technology. It is developed on a foundation of MC38 cells with a knockout of mouse PDL1.
<b>Quantity</b>	5E6 Cells per vial, 1 mL
<b>Product Format</b>	1 vials of frozen cells
<b>Shipping</b>	Shipped on dry ice
<b>Storage Conditions</b>	Liquid nitrogen immediately upon receipt
<b>Target</b>	Human_PDL1
<b>Gene ID/Uniprot ID</b>	NP_054862.1
<b>Host Cell</b>	MC38
<b>Recovery Medium</b>	DMEM+10% FBS+1% P.S
<b>Growth medium</b>	DMEM+10% FBS+1% P.S+2 µg/mL Blasticidin+200 µg/mL Hygromycin+2.5 µg/mL Puromycin
<b>Note</b>	None
<b>Freezing Medium</b>	90% FBS+10% DMSO
<b>Growth properties</b>	Adherent
<b>Growth Conditions</b>	37°C, 5% CO <sub>2</sub>
<b>Mycoplasma Testing</b>	The cell line has been screened to confirm the absence of Mycoplasma species.
<b>Safety considerations</b>	Biosafety Level 2
<b>Note</b>	It is recommended to expand the cell culture and store a minimum of 10 vials at an early passage for potential future use.

## Materials

Reagent	Manufacturer/Catalogue No.
DMEM	VivaCell/C3110-0500
Fetal Bovine Serum	ExCell/FSP500
Pen/Strep	Thermo/15140-122
Blasticidin	Genomeditech/GM-040404
Hygromycin	Genomeditech/GM-040403
Puromycin	Genomeditech/GM-040401
Anti-PDL1 hIgG4 Reference Antibody(Adebio)	Genomeditech/GM-86853MAB
PE anti-mouse CD274 (B7-H1, PD-L1) Antibody	BioLegend/124307

## Figures




Figure 1 | H\_PDL1 MC38(mouse PDL1 KO) Cell Line (Cat. GM-C20064) was determined by flow cytometry using Anti-PDL1 hIgG4 Reference Antibody (Adebio) (Cat. [GM-86853MAB](#)).





Figure 2 | H\_PDL1 MC38(mouse PDL1 KO) Cell Line (Cat. GM-C20064) was determined by flow cytometry using PE anti-mouse CD274 (B7-H1, PD-L1) Antibody(BioLegend/124307).






Figure 5 | Body weight changes after implantation of H\_PDL1 MC38(mouse PDL1 KO) in C57BL/6 mice. Under the same conditions, body weight was measured twice per week and is presented as mean  $\pm$  SEM.

## Cell Recovery

Recovery Medium: DMEM+10% FBS+1% P.S

To insure the highest level of viability, thaw the vial and initiate the culture as soon as possible upon receipt. If upon arrival, continued storage of the frozen culture is necessary, it should be stored in liquid nitrogen vapor phase and not at  $-70^{\circ}\text{C}$ . Storage at  $-70^{\circ}\text{C}$  will result in loss of viability.

- Thaw the vial by gentle agitation in a  $37^{\circ}\text{C}$  water bath. To reduce the possibility of contamination, keep the O-ring and cap out of the water. Thawing should be rapid (approximately 2 - 3 minutes).
- Remove the vial from the water bath as soon as the contents are thawed, and decontaminate by dipping in or spraying with 70% ethanol. All of the operations from this point on should be carried out under strict aseptic conditions.
- Transfer the vial contents to a centrifuge tube containing 5.0 mL complete culture medium and spin at approximately  $176 \times g$  for 5 minutes. Discard supernatant.
- Resuspend cell pellet with the recommended recovery medium. And dispense into appropriate culture dishes.
- Incubate the culture at  $37^{\circ}\text{C}$  in a suitable incubator. A 5%  $\text{CO}_2$  in air atmosphere is recommended if using the medium described on this product sheet.

## Cell Freezing

Freezing Medium: 90% FBS+10% DMSO

- Centrifuge at  $176 \times g$  for 3 minutes to collect cells.
- Resuspend the cells in pre-cooled freezing medium and adjust the cell density to  $5 \times 10^6$  cells/mL.
- Aliquot 1 mL into each vial.

- d) Place the vial in a controlled-rate freezing container and store at -80°C for at least 1 day, then transfer to liquid nitrogen as soon as possible.

## Cell passage

Growth medium: DMEM+10% FBS+1% P.S+2 µg/mL Blasticidin+200 µg/mL Hygromycin+2.5 µg/mL Puromycin

For the first 1 to 2 passages post-resuscitation, use the recovery medium. Once the cells have stabilized, switch to a growth medium.

- Remove and discard culture medium.
- Briefly rinse the cell layer with PBS to remove all traces of serum that contains trypsin inhibitor.
- Add 1.0 mL of 0.25% (w/v) Trypsin-EDTA solution to dish and observe cells under an inverted microscope until cell layer is dispersed (usually within 30 to 60 seconds at 37°C).
- Note: To avoid clumping do not agitate the cells by hitting or shaking the flask while waiting for the cells to detach. Cells that are difficult to detach may be placed at 37°C to facilitate dispersal.
- Add 2.0 mL of growth medium to mix well and aspirate cells by gently pipetting.
- After centrifugation, resuspend the pellet and add appropriate aliquots of the cell suspension to new culture vessels.
- Incubate cultures at 37°C.

**Subcultivation Ratio: A subcultivation ratio of 1:4 - 1:5 is recommended**

**Medium Renewal: Every 2 to 3 days**

## Notes

- After the stabilization of the cell condition, there will be fewer dead cells post-passage, the cell growth rate will tend to stabilize, cell morphology will become uniform, and the cells will appear robust.

## Sequence

CD274(PD-L1) [NP\\_054862.1](#)

MRIFAVFIFMTYWHLNNAFTVTVPKDLYVVEYGSNMTECKFPVEKQLDLAALIVYWEMEDKNIIQFVHGEE  
DLKVQHSSYRQRARLLKQQLSLGNAALQITDVKLQDAGVYRCMISYGGADYKRITVKVNAPYNKINQRILV  
VDPVTSEHELTCQAEGYPKAEVIWTSSDHQVLSGKTTTTNSKREEKLFNVTSTLRINTTTNEIFYCTFRRLDPE  
ENHTAELVIPELPLAHPNERTHLVILGAILLCLGVALTFIFRLRKGRMMDVKKCGIQDTSKKQSDTHLEET\*

## Related Products

PD-1:PD-L1(B7-H1):PDL2	
<a href="#">Mouse_PDL1 KO CT26 Cell Line</a>	<a href="#">Mouse_PDL1 KO LLC1 Cell Line</a>
<a href="#">Mouse_PDL1 KO MC38 Cell Line</a>	<a href="#">aAPC(OKT3) PDL1 CHO-K1 Cell Line</a>
<a href="#">H_PD-1 Reporter Jurkat Cell Line</a>	<a href="#">H_PD1 SHP1 Reporter Jurkat Cell line</a>
<a href="#">H_PD1 SHP2 Reporter Jurkat Cell Line</a>	<a href="#">H_PDCD1LG2(PDL2) aAPC CHO-K1 Cell Line</a>
<a href="#">Mouse PDL1 aAPC CHO-K1 Cell Line</a>	<a href="#">Mouse_PD-1 Reporter Jurkat Cell Line</a>

Canine_PD-1 CHO-K1 Cell Line	Canine_PD-1 HEK-293 Cell Line
Cynomolgus_PD1 CHO-K1 Cell Line	Cynomolgus_PD-L1 HEK-293 Cell Line
H_CD274(PD-L1) CHO-K1 Cell Line	H_CD274(PD-L1) MC38 Cell Line
H_PDCD1(PD-1) CHO-K1 Cell Line	H_PDCD1(PD-1) CHO-K1 Cell Line (Low Expression)
H_PDCD1(PD-1) HEK-293 Cell Line	H_PDCD1LG2(PDL2) CHO-K1 Cell Line
H_PDL1 CT26(mouse PDL1 KO) Cell Line	H_PD-L1 HEK-293 Cell Line
H_PDL1 LLC1(mouse_PDL1 KO) Cell Line	H_PDL1 LLC1(mouse_PDL1 KO) Cell Line
H_PD-L1 Raji Cell Line	M_PDCD1(PD-1) CHO-K1 Cell Line
Anti-Canine_PD1 mIgG2a Antibody(4F12-E6)	Anti-CTLA-4/PD-1 hIgG1 Bispecific Antibody(Cadonilimab)
Anti-CTLA4×PD-1 hIgG1 Reference Antibody (Cadbio)	Anti-H_CD274(PDL1) hIgG1 Antibody(Atezolizumab)
Anti-H_PDCD1(PD1) hIgG1 Antibody(Budigalimab)	Anti-H_PDCD1LG2 mIgG1 Antibody(3G2)
Anti-H_PDL1 hIgG1 Reference Antibody(Atezbio)	Anti-mouse PD1 RIGG2a Antibody(RMP1-14)
Anti-mouse PD-L1 mIgG1 Antibody(10F.9G2)	Anti-Mouse_PD1 mIgG1 Antibody(29F.1A12)
Anti-mouse_PD1 mIgG1 Antibody(RMP1-14)	Anti-Mouse_PD1×VEGF hIgG1 Bispecific Antibody
Anti-PD1 hIgG1 Antibody(Peresolimab)	Anti-PD1 hIgG1 Reference Antibody (Perbio)
Anti-PD1 hIgG1 Reference Antibody(Rosnbio)	Anti-PD1 hIgG4 Antibody(Pembrolizumab)
Anti-PD1 hIgG4 Reference Antibody (Nivbio)	Anti-PD1 hIgG4 Reference Antibody (Pembio)
Anti-PD1 hIgG4 Reference Antibody (Sintbio)	Anti-PD-1 hIgG4 Reference Antibody (Torbio)
Anti-PD1 hIgG4 Reference Antibody(Cambio)	Anti-PD-1 hIgG4 Reference Antibody(Tislbio)
Anti-PD1-IL2v Fusion hIgG1 Antibody(2149)	Anti-PD1-IL2v Fusion hIgG1 Antibody(KY-0118)
Anti-PD-L1 hIgG1 Reference Antibody(Avebio)	Anti-PDL1 hIgG4 Reference Antibody(Adebio)
Anti-PD-L2 hIgG1 Antibody(Hz25G4-1.1)	Anti-VEGF×PD1 hIgG1 Reference Antibody (Ivobio)
Biotinylated Human PD1 Protein; His-Avi Tag	Biotinylated Human PDL1 Protein; His-Avi Tag
Canine PD1 Protein; hFc Tag	Cynomolgus PDL1 Protein; His Tag
Human PD1 Protein; hFc Tag	Human PD1 Protein; His Tag
Human PDL1 Protein; His Tag	Human PDL1 Protein; mFc Tag
Human PDL2 Protein; mFc Tag	Mouse PDL1 Protein; His Tag
<b>In Vivo MAb Isotype Controls</b>	
Biotinylated Human IgE Isotype Control; His-Avi Tag (Anti-RSV)	Human IgG1 (K214R) Isotype Control(Anti-HEL)
Human IgG1 (K214R;D356E;L358M) Isotype Control(Anti-HEL)	Human IgG1 (KRLALAGA) Isotype Control(Anti-HEL)
Human IgG1 Isotype Control(Anti-HEL)	Human IgG1 Isotype Control(Anti-RSV)
Human IgG1 Isotype Control(MOPC-21)	Human IgG1(LALA) Isotype Control(Anti-HEL)
Human IgG1(LALAPG) Isotype Control(Anti-HEL)	Human IgG1(LALAPG;N297A) Isotype Control(Anti-HEL)
Human IgG1(N297A) Isotype Control(Anti-HEL)	Human IgG4(S228P) Isotype Control(Anti-HEL)
Mouse IgG1 (D265A) Isotype Control(Anti-HEL)	Mouse IgG1 Isotype Control(Anti-HEL)
Mouse IgG2a Isotype Control(Anti-HEL)	Mouse IgG2a Isotype Control(Anti-RSV)
Mouse IgG2a(D265A) Isotype Control(Anti-HEL)	Mouse IgG2b Isotype Control(Anti-HEL)
Rat IgG1 Isotype Control(Anti-HEL)	

## License Agreement:

**By purchasing and using this cell line product, the user voluntarily agrees to accept and abide by the following policies:**

- This cell line product is restricted to research use only and shall not be used for any commercial purposes.
- This product is strictly prohibited from being used in the diagnosis or treatment of human or animal diseases, and shall not be directly used in experiments involving humans.
- Users and their contractors engaged for their benefit may use this material and its derivatives only within the agreed research scope; modification of the material is not permitted, nor may it be distributed, sold, transferred, or otherwise provided to any other entity (including affiliates).
- If use beyond the above scope is required, prior written permission from Genomeditech (Shanghai) Co.,Ltd. must be obtained. For details, please contact Genomeditech (Shanghai) Co.,Ltd.

Genomeditech