


pET-41a(+)

pET41 系列载体是设计用来克隆和高水平表达蛋白质的载体，载体融合表达 GST 标签蛋白纯化片段。T7 蛋白表达区在质粒图谱上面是反向的。

图谱信息：



质粒元件信息：

f1 ori	12..467
KanR	560..1375
ori	1497..2085
bom	2271..2413
rop	2515..2706
lacI	3278..4360
lacI promoter	4361..4438
T7 promoter	4751..4769
lac operator	4770..4794
GST	4839..5492
6xHis	5520..5537
thrombin site	5547..5564
S-Tag	5580..5624
enterokinase site	5655..5669
8xHis	5761..5784
T7 terminator	5861..5908

