Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Anti-H_Nectin4 hlgG1 Antibody(Enfortumab)

Product information

GM-23369AB-10 10 μg GM-23369AB-100 100 μg GM-23369AB-1000 1 mg

Antibody Information

Species Reactivity Human; Cynomolgus

Clone Enfortumab

Source/Isotype Monoclonal human IgG1

Application Flow cytometry
Specificity Detects Nectin4

Gene Nectin4

Other Names EDSS1, LNIR, PRR4, PVRL4, nectin-4
Gene ID 81607(human), 102127936(cynomolgus)

Background This gene encodes a member of the nectin family. The encoded protein

contains two immunoglobulin-like (Ig-like) C2-type domains and one Ig-like V-type domain. It is involved in cell adhesion through trans-homophilic and -heterophilic interactions. It is a single-pass type I membrane protein. The soluble form is produced by proteolytic cleavage at the cell surface by the metalloproteinase ADAM17/TACE. The secreted form is found in both breast tumor cell lines and breast tumor patients. Mutations in this gene are the cause of ectodermal dysplasia-syndactyly syndrome type 1,

an autosomal recessive disorder.

Storage Store at 2-8°C short term (1-2 weeks). Store at ≤ -20°C long term. Avoid

repeated freeze-thaw.

Formulation Phosphate-buffered solution, pH 7.2.

Endotoxin < 1 EU/mg, determined by LAL gel clotting assay

Version:3.1 Revision Date:12/25/2023

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825


Toll-free: 400 627 9288

Email: service@genomeditech.com

Data Examples

Flow cytometry

The recommended usage range is 0.5-4 μ g per test. H_NECTIN4 HEK-293 Cell Line (Catalog # GM-C19280)was stained with Anti-H_Nectin4 hlgG1 Antibody (Catalog # GM-23369AB) or isotype control antibody, followed by anti-Human IgG APC-conjugated Secondary Antibody .




L	SampleID	Geometric Mean : RL1-H
I	HEK-293 anti-NECTIN4+APC-2nd Ab	152
	HEK-293 H_NECTIN4 H_IgG+APC-2nd Ab	164
I	HEK-293 H_NECTIN4 anti-NECTIN4+APC-2nd Ab	86016

Fig. FACS

Flow cytometry

The recommended usage range is 0.5-4 μ g per test. Cynomolgus_Nectin4 CHO-K1 Cell Line(Catalog # GM-C19258) was stained with Anti-H_Nectin4 hlgG1 Antibody (Catalog # GM-23369AB) or isotype control antibody, followed by anti-Human lgG APC-conjugated Secondary Antibody .



L	SampleID	Geometric Mean : RL1-H
I	CHO-K1 anti-H_Nectin4+APC-2nd Ab	112
l	CHO-K1 Cyno_Nectin4 H_lgG+APC-2nd Ab	99.6
Į	CHO-K1 Cyno_Nectin4 anti-H_Nectin4+APC-2nd Ab	33596

Fig. FACS