

Anti-H_HLA-G1 hlgG1 Antibody(38373)

Product information

GM-28208AB-10	10 µg
GM-28208AB-100	100 µg
GM-28208AB-1000	1 mg

Antibody Information

Species Reactivity	Human
Clone	38373
Source/Isotype	Monoclonal human IgG1, κ
Application	Flow cytometry
Other Names	HLA-G, MHC-G
Gene ID	P17693-1 (human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Flow cytometry H_HLA-G1 CHO-K1 Cell Line (Catalog # GM-C15840) was stained with Anti-H_HLA-G1 hlgG1 Antibody(38373) or isotype control antibody.




Fig. FACS

Version:0.3.2.240521