

## Anti-ACE2 hlgG4 Antibody(3E8)

### Product information


GM-32019-10	10 µg
GM-32019-100	100 µg
GM-32019-1000	1 mg

### Antibody Information

Species Reactivity	Human
Clone	3E8
Source/Isotype	Monoclonal human IgG4, κ
Application	Flow cytometry
Other Names	ACEH
Gene ID	59272 (human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

### Data Examples

Flow cytometry H\_ACE2 HEK-293T Cell Line (Catalog # GM-C12968)was stained with Anti-ACE2 hlgG4 Antibody(3E8) or isotype control antibody.



SampleID	Geometric Mean : FL11-H
HEK-293T anti-ACE2+APC-2nd Ab	1954
HEK-293T H_ACE2 H_IgG+APC-2nd Ab	1193
HEK-293T H_ACE2 anti-ACE2+APC-2nd Ab	5.17E5

Fig. FACS

Version:0.3.2.240521