

Anti-H_KLRC1 hlgG4 Antibody(Monalizumab)

Product information


GM-34467AB-10	10 µg
GM-34467AB-100	100 µg
GM-34467AB-1000	1 mg

Antibody Information

Species Reactivity	Human
Clone	Monalizumab
Source/Isotype	Monoclonal human IgG4, κ
Application	Flow cytometry;
Other Names	CD159A, NKG2, NKG2A
Gene ID	NP_998823.2 (human);
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Flow cytometry H_NKG2A H_CD94 CHO-K1 Cell Line (Catalog # GM-C15621) was stained with Anti-H_KLRC1 hlgG4 Antibody(Monalizumab) or isotype control antibody.



SampleID	Geometric Mean : FL11-H
CHO-K1 anti-H_KLRC1(NKG2A)+APC-2nd Ab	1571
CHO-K1 H_NKG2A-CD94 H_IgG+APC-2nd Ab	947
CHO-K1 H_NKG2A-CD94 anti-H_KLRC1(NKG2A)+APC-2nd Ab	138423

Fig. FACS

Version:0.3.2.240521