

Anti-H_IL-22R1 hlgG1 Antibody(ARGX-112)

Product information

GM-46382AB-10	10 µg
GM-46382AB-100	100 µg
GM-46382AB-1000	1 mg

Antibody Information

Species Reactivity	Human
Clone	ARGX-112
Source/Isotype	Monoclonal human IgG1, λ
Application	Flow cytometry
Other Names	IL-21, IL-22, IL-D110, IL-TIF, ILTIF, TIFIL-23, TIFa, zcyto18
Gene ID	Q8N6P7-1(human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Flow cytometry H_IL22 Reporter 293 Cell Line(Catalog # GM-C19682) was stained with Anti-H_IL-22R1 hlgG1 Antibody(ARGX-112) or isotype control antibody.




Fig. FACS

Version:0.3.2.240521