

Anti-H_IL-22 hlgG1 Antibody(Fezakinumab)

Product information

GM-46509AB-10	10 µg
GM-46509AB-100	100 µg
GM-46509AB-1000	1 mg

Antibody Information

Species Reactivity	Human
Clone	Fezakinumab
Source/Isotype	Monoclonal human IgG1,λ
Application	Block assay
Other Names	IL-21, IL-22, IL-D110, IL-TIF, ILTIF, TIFIL-23, TIFa, zcyto18
Gene ID	Q8N6P7-1(human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Block assay Anti-H_IL-22 hlgG1 Antibody inhibits H_IL22 Reporter 293 Cell Line (Catalog # GM-C19682) Luminescence induced by H_IL22 Protein.




Fig. assay