

## Anti-H\_CB1 hlgG4 Antibody(Nimacimab)

### Product information

GM-47472AB-10	10 µg
GM-47472AB-100	100 µg
GM-47472AB-1000	1 mg

### Antibody Information

Species Reactivity	Human
Clone	Nimacimab
Source/Isotype	Monoclonal human IgG4, κ
Application	Flow cytometry;
Other Names	CANN6, CB-R, CNR1, CB1A, CB1K5, CB1R, CNR
Gene ID	P21554-1 (human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

### Data Examples

Flow cytometry H\_CB1 HEK-293 Cell Line (Catalog # GM-C22609) was stained with Anti-H\_CB1 hlgG4 Antibody(Nimacimab) or isotype control antibody.




Fig. FACS

Version:0.3.2.240521