

# Anti-OSMR hIgG4 Antibody(Vixarelimab)

## Product information

GM-50874AB-10	10 µg
GM-50874AB-100	100 µg
GM-50874AB-1000	1 mg

## Antibody Information

Species Reactivity	Human
Clone	Vixarelimab
Source/Isotype	Monoclonal human hIgG4, λ
Application	Flow cytometry;
Other Names	OSMRB; PLCA1; IL-31RB; OSMRbeta; IL-31R-beta
Gene ID	Q8NI17-1(human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

## Data Examples

Flow cytometry H\_IL-31 Reporter Cell Line (Catalog # GM-C25989) was stained with Anti-OSMR hIgG4 Antibody(Vixarelimab) or isotype control antibody.




Fig. FACS