

Anti-IL-1 β hlgG1 Antibody(Canakinumab)

Product information

GM-51764AB-10	10 μ g
GM-51764AB-100	100 μ g
GM-51764AB-1000	1 mg

Antibody Information

Species Reactivity	Human;
Clone	Canakinumab
Source/Isotype	Monoclonal human IgG1, κ
Application	Block assay;
Other Names	IL-1, IL1-BETA, IL1F2, IL1beta, IL1B
Gene ID	3553(human)
Storage	Store at 2-8 $^{\circ}$ C short term (1-2 weeks).Store at \leq -20 $^{\circ}$ C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Block assay Anti-IL-1 β hlgG1 Antibody(Canakinumab) inhibits H_IL-1 β Reporter 293 Cell Line (Catalog # GM-C26463) Luminescence induced by Human IL-1 beta / IL1B Protein.




Fig. assay

Version:0.3.2.240521