

Anti-CD132(IL2RG) hIgG4 Antibody(REGN7257)

Product information

GM-52334AB-10	10 µg
GM-52334AB-100	100 µg
GM-52334AB-1000	1 mg

Antibody Information

Species Reactivity	Human;
Clone	REGN7257
Source/Isotype	Monoclonal human IgG4, κ
Application	Flow cytometry;
Other Names	CIDX; IL-2RG; IMD4; P64; SCIDX; SCIDX1
Gene ID	NM_000206.1 (human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Flow cytometry

H_CD122 CD132 Reporter Cell Line (Catalog # GM-C29054) was stained with Anti-CD132(IL2RG) hIgG4 Antibody(REGN7257) or isotype control antibody.




Fig. FACS