

## Anti-MICA/MICB mIgG2a Antibody(PDI-01)

### Product information


GM-52540AB-10	10 µg
GM-52540AB-100	100 µg
GM-52540AB-1000	1 mg

### Antibody Information

Species Reactivity	Cynomolgus
Clone	PDI-01
Source/Isotype	Monoclonal Mouse IgG2a, κ
Application	Flow cytometry;
Other Names	MIC-A; PERB11.1
Gene ID	Q2MGE0-1(cynomolgus)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

### Data Examples

Flow cytometry Cynomolgus\_MICA(Q2MGE0-1) CHO-K1 Cell Line (Catalog # GM-C26177) was stained with Anti-MICA/MICB mIgG2a Antibody(PDI-01) or isotype control antibody.



SampleID	Geometric Mean : FL11-H
CHO-K1 anti-MICA/MICB+APC-2nd Ab	717
CHO-K1 Cyno_MICA M_IgG+APC-2nd Ab	1233
CHO-K1 Cyno_MICA anti-MICA/MICB+APC-2nd Ab	207689

Fig. FACS

Version:0.3.2.240521