Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Mouse CDH17 Protein; His Tag

Product Information

Product Name Mouse CDH17 Protein; His Tag

Storage temp. Store at \leq -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-75418RP-100 / 100 μg

GM-75418RP-1000 / 1 mg

Protein Information

Alternative Names Cadherin-17,HPT-1,LI-cadherin

Source Mouse CDH17 Protein; His Tag (GM-75418RP) is expressed from human 293

cells (HEK-293). It contains AA Phe26-Met786 (Accession # Q9R100-1).

This protein carries a His tag at the C-terminus.

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/μg, determined by LAL gel clotting assay

Predicted Mol Mass 86.3 KDa

Formulation Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

Description Cadherin-17 is a protein that in humans is encoded by the CDH17 gene. This gene

is a member of the cadherin superfamily, genes encoding calcium-dependent, membrane-associated glycoproteins. The encoded protein is cadherin-like, consisting of an extracellular region, containing 7 cadherin domains, and a transmembrane region but lacking the conserved cytoplasmic domain. The protein is a component of the gastrointestinal tract and pancreatic ducts, acting as an intestinal proton-dependent peptide transporter in the first step in oral absorption of many medically important peptide-based drugs. The protein may also play a role in the morphological organization of liver and intestine.


Version:3.3 Revision Date:25/12/2023

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


SDS-PAGE



This protein on SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

Bioactivity-ELISA

Mouse CDH17 Protein; His Tag, ELISA 3 μg of Monoclonal Anti-CDH17 Antibody per well



Anti-CDH17 Antibody (Catalog # GM-75321AB) was immobilized at 30 μ g/ml (100 μ L/well). Increasing concentrations of Mouse CDH17 Protein; His Tag (Catalog # GM-75418RP) were added.