Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

# Human IgG4(S228P) Isotype Control(Anti-HEL)

#### **Product Information**

Product Name Human IgG4(S228P) Isotype Control(Anti-HEL)

**Storage temp.** Store at  $2-8^{\circ}$ C short term (1-2 weeks). Store at  $\leq -20^{\circ}$ C long term. Avoid

repeated freeze-thaw.

Catalog# / Size GM-77406AB-1mg / 1 mg

GM-77406AB-5mg / 5 mg GM-77406AB-25mg / 25 mg GM-77406AB-50mg / 50 mg GM-77406AB-100mg / 100 mg

#### **Antibody Information**

Expression System CHO

Aggregation < 5% as determined by SEC-HPLC

Purity > 95% as determined by SDS-PAGE

Endotoxin <1 EU/mg, determined by LAL gel clotting assay

Sterility 0.22 µm Filtered

Target HEL Clone /

Alternative Names /

Source/Isotype Monoclonal Human IgG4 S228P, kappa

Application Negative control

**Description**A homologous control antibody is an antibody that is used to control the

target antibody in an experiment. The homologous control antibody should have the same immunoglobulin type and concentration as the test antibody, but not have the same specificity as the target antibody. It is often used to help researchers distinguish between non-specific signals and background signals in experiments to ensure the accuracy and

reliability of experimental results.


Formulation PBS, pH 7.4, Contains no stabilizers or preservatives

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288


Email: service@genomeditech.com

### **SDS-PAGE**



On SDS-PAGE under reducing (R)/non-reducing(N-R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

## **SEC-HPLC**



The purity of this product is more than 95% verified by SEC-HPLC.