Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

## Human OX40L Protein; mFc Tag

#### **Product Information**

Product Name Human OX40L Protein; mFc Tag

Storage temp. Store at  $\leq$  -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-83969RP-100 / 100 μg

GM-83969RP-1000 / 1 mg

#### **Protein Information**

Alternative Names TNFSF4,CD252,Glycoprotein Gp34,TXGP1,CD134 ligand,CD134L

Source Human OX40L Protein; mFc Tag (GM-83969RP) is expressed from human 293

cells (HEK-293). It contains AA Gln51-Leu183 (Accession # P23510-1).

This protein carries a mFc tag at the C-terminus.

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/μg, determined by LAL gel clotting assay

Predicted Mol Mass 42.3 KDa

**Formulation** Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

**Description** OX40L is a cell surface molecule, also known as CD252. It is a receptor ligand

and a member of the tumor necrosis factor family. OX40L typically binds to its receptor OX40, thereby transmitting signals between immune cells and regulating the activation and proliferation of T cells. This signaling is crucial for the normal functioning of the immune system, making OX40L potentially important in


immune regulation and therapy.

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


#### **SDS-PAGE**



On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

### **Bioactivity-ELISA**

# Human OX40L Protein; mFc Tag, ELISA 0.2 μg of Human OX40L Protein; mFc Tag per well



Human OX40L Protein; mFc Tag (Catalog # GM-83969RP) was immobilized at 2  $\mu$ g/ml (100  $\mu$ L/well). Increasing concentrations of Monoclonal Anti-OX40L hIgG4 Antibody (Amlitelimab) (Catalog # GM-82533AB) were added.