Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

## Cynomolgus TL1A Protein; His Tag

#### **Product Information**

Product Name Cynomolgus TL1A Protein; His Tag

**Storage temp.** Store at  $\leq$  -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-84080RP-100 / 100 μg

GM-84080RP-1000 / 1 mg

#### **Protein Information**

Alternative Names VEGI, TNFSF15

Source Cynomolgus TL1A Protein; His Tag (GM-84080RP) is expressed from human

293 cells (HEK-293). It contains AA Leu 72 - Leu 251 (Accession # F6S8F9-1).

This protein carries a His tag at the N-terminus.

Purity > 90% as determined by SDS-PAGE

Endotoxin < 1 EU/μg, determined by LAL gel clotting assay

Predicted Mol Mass 22.3 KDa

**Formulation** Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

**Description** Tumor necrosis factor (TNF)-like cytokine 1A (TL1A) is a member of the TNF

superfamily of proteins (TNFSF15), which signals through association with death domain receptor 3 (DR3). Decoy receptor 3 (DcR3) competes with DR3 for TL1A binding and inhibits functional signaling. These proteins are significantly upregulated in inflamed intestinal tissues, and their pathogenetic importance for inflammatory bowel disease (IBD) is suggested by accumulating evidence. TL1A/DR3 induce costimulatory signals to activated lymphocytes, including the gut-specific populations of CD4+CD161+ and CD4+CCR9+ cells, affecting all major effector pathways and inducing the mucosal upregulation of Th1, Th2, and


Th17 factors.

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


#### **SDS-PAGE**



On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 90%.

#### **Bioactivity-ELISA**

# Cynomolgus TL1A Protein; His Tag, ELISA 0.5 µg of Cynomolgus TL1A Protein; His Tag per well



Anti-H\_TNFSF15(TL1A) hIgG1 Antibody(Tulisokibart、PRA-023) Conc.(ng/mL)

Cynomolgus TL1A Protein; His Tag (Catalog # GM-84080RP) was immobilized at 5  $\mu$ g/ml (100  $\mu$ L/well). Increasing concentrations of Monoclonal Anti-H\_TNFSF15(TL1A) hIgG1 Antibody(Tulisokibart、PRA-023) (Catalog # GM-58915AB) were added.


Version:3.3 Revision Date:25/12/2023

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

### **Bioactivity CELL BASE**



Cynomolgus TL1A Protein; His Tag (Catalog # GM-84080RP) was added to  $H_{TNFSF15}$ (TL1A) Reporter Cell Line (Catalog # GM-C30289), and stimulates TL1A/DR3 signals.