Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

## Cynomolgus B7H4 Protein; His Tag

#### **Product Information**

**Product Name** Cynomolgus B7H4 Protein; His Tag

**Storage temp.** Store at  $\leq$  -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-84641RP-100 / 100 μg

GM-84641RP-1000 / 1 mg

#### **Protein Information**

Alternative Names B7-H4,VTCN1,B7S1,B7h.5

Source Cynomolgus B7H4 Protein; His Tag (GM-84641RP) is expressed from human

293 cells (HEK-293). It contains AA Phe29-Ala258 (Accession # F7B770-1).

This protein carries a His tag at the C-terminus.

Purity > 90% as determined by SDS-PAGE

Endotoxin <1 EU/μg, determined by LAL gel clotting assay

Predicted Mol Mass 27.2 KDa

**Formulation** Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

Description B7H4 protein is a cell surface protein belonging to the B7/CD28 immune

regulatory family. This protein plays a crucial role in regulating immune responses and tumor immune evasion. Overexpression of B7H4 is associated with the development and progression of various tumors, making it a potential target for cancer therapy. Additionally, B7H4 may also promote tumor immune evasion by influencing the function of other immune cells in the tumor microenvironment, such as inhibiting natural killer cells (NK cells) and regulatory T cells (Treg cells). The B7H4 protein is of great interest in the fields of immunotherapy and tumor


immunotherapy.

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


#### **SDS-PAGE**



On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 90%.

### **Bioactivity-ELISA**

# Cynomolgus B7H4 Protein; His Tag ELISA 0.2 µg Cynomolgus B7H4 Protein; His Tag of per well



Cynomolgus B7H4 Protein; His Tag (Catalog # GM-84641RP) was immobilized at 2  $\mu$ g/ml (100  $\mu$ L/well). Increasing concentrations of Anti-H\_B7H4(VTCN1) hIgG1 Antibody (FPA150) (Catalog # GM-24028AB) were added.