Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Biotinylated Human ALPP Protein; His-Avi Tag

Product Information

Product Name Biotinylated Human ALPP Protein; His-Avi Tag

Storage temp. Store at \leq -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-85159RP-25 / 25 μg

GM-85159RP-200 / 200 µg

Protein Information

Alternative Names PLAP

Source Biotinylated Human ALPP Protein; His-Avi Tag (GM-85159RP) is expressed

from human 293 cells (HEK-293). It contains AA Ile23-Asp506 (Accession #

P05187-1).

This protein carries a His-Avi tag at the C-terminus.

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/μg, determined by LAL gel clotting assay

Predicted Mol Mass 57.4 KDa

Formulation Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

Description PLAP is the abbreviation for Placental alkaline phosphatase, which is an enzyme

protein. PLAP is highly expressed in placental tissue, hence its name. The human placental ALP (PLAP) gene was identified on chromosome 2 and possesses 87% similarity with IAP gene. Interestingly, cancer cells can re-express PLAP gene as Regan isoenzyme. In contrast to the other ALP isoenzymes, PLAP is a polymorphic enzyme with up to 18 allelozymes originating from point mutations. As a result of the low catalytic activity in comparison with the other isoenzymes, PLAP is considered to be inefficient due to the neutral pH of placenta. Furthermore, the structure of PLAP consists of two disulfide bridges: the first one stabilises the PLAP orientation, which is situated near the anchor of the Asp481 Glycosylyphosphatidylinositol (GPI) and the second stiffens the anchor of the carbohydrate chain and is situated close the first part of the glycosylation. There are several roles by PLAP, including the transmission of mother-to-fetus immunoglobulins G (IgG). PLAP stimulates DNA synthesis and fibroblast cell proliferation with insulin in the presence of Zn2+ and Ca2+. It is a synthesiser for a single tissue called syncytiotrophoblast and responsible for transfer of substances such as nutrients and oxygen. PLAP is an essential modulator of foetal

Version:3.3 Revision Date:25/12/2023


growth. PLAP is also a potential marker for various cancers.

Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


SDS-PAGE



On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

Bioactivity-ELISA

Biotinylated Human ALPP Protein; His-Avi Tag, ELISA 0.1 µg Biotinylated Human ALPP Protein; His-Avi Tag of per well



Biotinylated Human ALPP Protein; His-Avi Tag (Catalog # GM-85159RP) was immobilized at 1 μ g/ml (100 μ L/well) on streptavidin precoated. Increasing concentrations of Anti-ALPPL2 mIgG2a antibody (SGN-ALPV) (Catalog # GM-58773AB) were added.