

Biotinylated Human PDL1 Protein; His-Avi Tag

Product Information

Product Name	Biotinylated Human PDL1 Protein; His-Avi Tag
Storage temp.	Store at \leq -70°C, stable for 6 months after receipt. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.
Catalog# / Size	GM-85610RP-25 / 25 μg GM-85610RP-200 / 200 μg


Protein Information

Alternative Names	PD-L1,CD274,B7-H1,PDCD1L1,PDCD1LG1
Source	Biotinylated Human PDL1 Protein; His-Avi Tag (GM-85610RP) is expressed
	from human 293 cells (HEK-293). It contains AA Phe 19 - Arg 238 (Accession #
	NP_054862.1).
	This protein carries a His-Avi tag at the C-terminus.
Purity	> 95% as determined by SDS-PAGE
Endotoxin	< 1 EU/µg, determined by LAL gel clotting assay
Predicted Mol Mass	26.6 KDa
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH7.4.
Description	PDL1, also known as programmed death-ligand 1, is a protein that is expressed
	on the surface of certain cells, including cancer cells and immune cells. It plays a
	key role in regulating the immune response by binding to its receptor, PD1, on T
	cells and inhibiting their activity. This interaction can suppress the immune
	system's ability to recognize and attack cancer cells, allowing tumors to evade
	destruction.
	PDL1 has become an important target for cancer immunotherapy, as blocking the
	interaction between PDL1 and PD1 can help to unleash the immune system's
	ability to recognize and destroy cancer cells. Drugs that target PDL1 or PD1,
	known as checkpoint inhibitors, have shown promising results in treating various
	types of cancer, including melanoma, lung cancer, and bladder cancer.
	The expression of PDL1 in tumor cells is also used as a biomarker to predict the
	response to checkpoint inhibitor therapy. Patients with high levels of PDL1
	expression in their tumors may be more likely to benefit from treatment with
	checkpoint inhibitors. As a result, PDL1 testing has become an important part of
	the decision-making process for cancer treatment.

Version:3.3 Revision Date:25/12/2023




SDS-PAGE



On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

Biotinylated Human PDL1 Protein; His-Avi Tag, ELISA

Bioactivity-ELISA



Biotinylated Human PDL1 Protein; His-Avi Tag (Catalog # GM-85610RP) was immobilized at 1 μg/ml (100 μL/well) on streptavidin precoated. Increasing concentrations of Anti-H_CD274(PDL1) hIgG1 Antibody (Atezolizumab) (Catalog # GM-31740AB) were added.