Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Human ACVR2A Protein; His Tag

Product Information

Product Name Human ACVR2A Protein; His Tag

Storage temp. Store at \leq -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-87256RP-100 / 100 μg

GM-87256RP-1000 / 1 mg

Protein Information

Alternative Names Activin receptor type IIA, ACTRIIA

Source Human ACVR2A Protein; His Tag (GM-87256RP) is expressed from human 293

cells (HEK-293). It contains AA Ala20-Pro134 (Accession # P27037).

This protein carries a His tag at the C-terminus.

Purity > 95% as determined by SDS-PAGE

Endotoxin $< 1 \text{ EU/}\mu\text{g}$, determined by LAL gel clotting assay

Predicted Mol Mass 13.4 KDa

Formulation Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

Description The ACVR2A protein is a receptor-type serine/threonine kinase and belongs to

the bone morphogenetic protein receptor family. The ACVR2A protein plays an important role in cell signaling and regulates biological processes such as cell

proliferation, differentiation, and apoptosis.

The ACVR2A protein is a transmembrane protein that contains an extracellular domain, a transmembrane domain, and an intracellular kinase domain. During signal transduction, the ACVR2A protein binds to its ligand, activates its kinase activity, phosphorylates downstream signaling molecules, and triggers a series of signaling pathways that affect cell function.

The ACVR2A protein plays an important role in embryonic development, organ formation, bone growth and repair, immune regulation, and other physiological processes. It is also associated with the occurrence and development of various diseases such as cancer, inflammatory diseases, and autoimmune diseases. Therefore, the ACVR2A protein is considered an important therapeutic target and


can be used to develop drugs for related diseases.

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

SDS-PAGE




On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

Bioactivity-ELISA

Bioactivity-ELISA

0.1 μg Human ACVR2A Protein; His Tag of per well



Human ACVR2A Protein; His Tag (Catalog # GM-87256RP) was immobilized at 1 μ g/ml (100 μ L/well). Increasing concentrations of Anti-ACVR2B hIgG1 Antibody(Bimagrumab) (Catalog # GM-51148AB) were added.