Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Biotinylated Mouse ACVR2B Protein; His-Avi

Tag

Product Information

Product Name Biotinylated Mouse ACVR2B Protein; His-Avi Tag

Storage temp. Store at \leq -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-87261RP-25 / 25 μg

GM-87261RP-200 / 200 µg

Protein Information

Alternative Names ACVR2B, Activin RIIB, ACTRIIB, MGC116908

Source Biotinylated Mouse ACVR2B Protein; His-Avi Tag (GM-87261RP) is expressed

from human 293 cells (HEK-293). It contains AA Ser19-Thr134 (Accession #

NP 031423.1).

This protein carries a His-Avi tag at the C-terminus.

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/μg, determined by LAL gel clotting assay

Predicted Mol Mass 15.5 KDa

Formulation Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

Description Activin receptor type-2B (ACVR2B) is a transmembrane receptor protein that

belongs to the TGF-beta superfamily and is involved in the signaling pathways that regulate various cellular processes. It functions as a receptor for activin, a protein that plays a role in the regulation of cell growth, differentiation, and apoptosis. ACVR2B is particularly important in the context of skeletal and muscle development, as it is involved in the regulation of bone and muscle growth and repair. In addition to its role in normal physiological processes, ACVR2B has garnered significant interest in the field of regenerative medicine and muscle-related conditions. Research has shown that the manipulation of ACVR2B signaling pathways may have therapeutic potential for conditions such as muscle wasting, muscle degeneration, and muscle injuries. This has led to investigations into the development of potential treatments targeting ACVR2B for muscle-related disorders. Overall, ACVR2B is a protein with diverse roles in cellular signaling, development, and disease, and ongoing research continues to uncover

its potential as a therapeutic target for various conditions, including muscle-

related disorders and cancer.


Version:3.3 Revision Date:25/12/2023

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


SDS-PAGE



On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

Bioactivity-ELISA

Biotinylated Mouse ACVR2B Protein; His-Avi Tag, ELISA 0.5 µg Biotinylated Mouse ACVR2B Protein; His-Avi Tag of per well



Biotinylated Mouse ACVR2B Protein; His-Avi Tag (Catalog # GM-87261RP) was immobilized at 5 μ g/ml (100 μ L/well) on streptavidin precoated. Increasing concentrations of Anti-ACVR2B hIgG1 Antibody (Bimagrumab) (Catalog # GM-51148AB) were added.