Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Mouse CTLA4 Protein; His Tag

Product Information

Product Name Mouse CTLA4 Protein; His Tag

Storage temp. Store at \leq -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-87557RP-100 / 100 μg

GM-87557RP-1000 / 1 mg

Protein Information

Alternative Names CTLA-4,CD152

Source Mouse CTLA4 Protein; His Tag (GM-87557RP) is expressed from human 293

cells (HEK-293). It contains AA Glu36-Phe162 (Accession # NP_033973).

This protein carries a His tag at the C-terminus.

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/μg, determined by LAL gel clotting assay

Predicted Mol Mass 14.7 KDa

Formulation Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

Description CTLA-4 (Cytotoxic T-lymphocyte-associated protein 4) is a protein receptor that

plays a crucial role as an immune checkpoint in regulating immune responses. It is expressed on the surface of T cells and helps regulate the immune system to prevent it from attacking healthy cells. CTLA-4 is essential for maintaining

immune homeostasis and preventing autoimmune diseases.

CTLA-4 inhibits T cell activation by competitively binding to B7 molecules (CD80 and CD86). It is rapidly expressed in the early stages of T cell activation and inhibits T cell activation and proliferation by suppressing co-stimulatory signals. This inhibitory function helps to prevent excessive immune responses and

the development of autoimmune diseases.

CTLA-4 is of significant importance in immunotherapy as inhibiting CTLA-4 can enhance immune responses and improve anti-tumor immune responses. Anti-CTLA-4 antibodies have been used in the treatment of various tumors such as melanoma. A deep understanding of the CTLA-4 protein and its function aids in the development of new immunotherapy strategies to enhance treatment

outcomes for diseases like cancer.


Version:3.3 Revision Date:25/12/2023

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


SDS-PAGE



On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

Bioactivity-ELISA

Mouse CTLA4 Protein; His Tag, ELISA 0.1 µg Mouse CTLA4 Protein; His Tag of per well



Mouse CTLA4 Protein; His Tag (Catalog # GM-87557RP) was immobilized at 1 μ g/ml (100 μ L/well). Increasing concentrations of Anti-CTLA4 mIgG2a Reference Antibody(9D9) (Catalog # GM-87422MAB) were added.