

Human EGFR Protein; His Tag

Product Information


Product Name	Human EGFR Protein; His Tag
Storage temp.	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.
Catalog# / Size	GM-87561RP-100 / 100 μg GM-87561RP-1000 / 1 mg

Protein Information

Alternative Names	EGF R,ERBB,ERBB1,HER1,PIG61
Source	Human EGFR Protein; His Tag (GM-87561RP) is expressed from human 293 cells (HEK-293). It contains AA Leu 25 - Ser 645 (Accession # P00533-1). This protein carries a His tag at the C-terminus.
Purity	> 95% as determined by SDS-PAGE
Endotoxin	< 1 EU/ μg , determined by LAL gel clotting assay
Predicted Mol Mass	69.0 KDa
Formulation	Supplied as a 0.2 μm filtered solution of PBS, pH7.4.
Description	<p>EGFR (Epidermal Growth Factor Receptor) protein is a receptor tyrosine kinase (RTK), belonging to the tyrosine kinase receptor family. EGFR is widely expressed in various cell types, including epithelial cells, nerve cells, and muscle cells. Activation of EGFR typically occurs through binding to its ligand (such as EGF), triggering its autophosphorylation and activating its tyrosine kinase activity. Activated EGFR then initiates a series of signaling pathways, including Ras-MAPK, PI3K-Akt, among others, thereby influencing cell physiological functions.</p> <p>EGFR plays a significant role in many types of cancers because abnormal EGFR signaling pathways can lead to excessive cell proliferation and malignant transformation. Therefore, EGFR has become one of the important targets for cancer treatment, with anti-cancer drugs targeting EGFR, such as EGFR inhibitors, being widely researched and applied in clinical therapy. EGFR holds importance in cell biology and cancer research, playing a crucial role in understanding cell signaling and cancer development mechanisms.</p>

Version:3.3 Revision Date:25/12/2023


SDS-PAGE



On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

Bioactivity-ELISA

Human EGFR Protein; His Tag, ELISA 0.1 μ g Human EGFR Protein; His Tag of per well



Human EGFR Protein; His Tag (Catalog # GM-87561RP) was immobilized at 1 μ g/ml (100 μ L/well). Increasing concentrations of Anti-H_EGFR hIgG1 Antibody(Necitumumab) (Catalog # GM-49155AB) were added.