Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Biotinylated Human HER3 Protein; His-Avi Tag

Product Information

Product NameBiotinylated Human HER3 Protein; His-Avi Tag

Storage temp. Store at \leq -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-87576RP-25 / 25 μg

GM-87576RP-200 / 200 µg

Protein Information

Alternative Names ERBB3,HER3,LCCS2,MDA-BF-1,MGC88033,c-erbB3,erbB3-S,p180-ErbB3

Source Biotinylated Human HER3 Protein; His-Avi Tag (GM-87576RP) is expressed

from human 293 cells (HEK-293). It contains AA Ser 20 - Thr 643 (Accession #

P21860-1)

This protein carries a His-Avi tag at the C-terminus.

Purity > 90% as determined by SDS-PAGE

Endotoxin < 1 EU/μg, determined by LAL gel clotting assay

Predicted Mol Mass 71.1 KDa

Formulation Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

Description HER3 protein, the full name of which is Human Epidermal Growth Factor

Receptor 3, is a member of the receptor tyrosine kinase (RTK) family and is involved in biological processes such as cell growth, proliferation, and differentiation. The HER3 protein is encoded by the ERBB3 gene and is part of the EGFR family, which includes EGFR (HER1), HER2, and HER4. Unlike other members of the EGFR family, HER3 protein has weak kinase activity itself and usually needs to form a dimer with other members to activate downstream

signaling pathways.

HER3 protein plays an important role in various cancers, including breast cancer, ovarian cancer, gastric cancer, and colorectal cancer. It is considered an important target for cancer treatment because overexpression or abnormal activation of HER3 is associated with tumor development and treatment resistance. Antibodies and inhibitors targeting HER3 have been developed and have shown some

therapeutic effects in clinical trials.

Overall, HER3 protein is an important cell signaling molecule with potential


importance in cancer research and treatment.

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

SDS-PAGE




On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 90%.

Bioactivity-ELISA

Bioactivity-ELISA

0.1 µg Biotinylated Human HER3 Protein; His-Avi Tag of per well



Anti-H ERBB3(HER3) hIgG1 Antibody(Barecetamab).(ng/mL)

Biotinylated Human HER3 Protein; His-Avi Tag (Catalog # GM-87576RP) was immobilized at 1 μ g/ml (100 μ L/well) on streptavidin precoated. Increasing concentrations of Anti-H_ERBB3(HER3) hIgG1 Antibody (Barecetamab) (Catalog # GM-28858AB) were added.