

Biotinylated Cynomolgus ITPRIPL1 Protein; His-Avi Tag

Product Information

Product Name Biotinylated Cynomolgus ITPRIPL1 Protein; His-Avi Tag

Storage temp. Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.
Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size **GM-87669RP-25 / 25 μg**
GM-87669RP-200 / 200 μg

Protein Information

Alternative Names ITPRIP like 1; ITPRIPL1; KIAA1754L; CD3L1

Source Biotinylated Cynomolgus ITPRIPL1 Protein; His-Avi Tag (GM-87669RP) is expressed from human 293 cells (HEK-293). It contains AA His 25 - Gly 103 (Accession # XP_005575041.2).

This protein carries a His-Avi tag at the C-terminus.

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/ μg , determined by LAL gel clotting assay

Predicted Mol Mass 11.2 KDa

Formulation Supplied as a 0.2 μm filtered solution of PBS, pH7.4.


Description ITPRIP like 1 (Inositol 1,4,5-trisphosphate receptor interacting protein-like 1) protein is a protein related to the inositol trisphosphate receptor, which is associated with the release of calcium ions stored within cells. This protein may play a role in regulating the release of calcium ions within cells and has an important role in cell signaling.

ITPRIP like 1 protein interacts with the inositol trisphosphate receptor (IP3R) to regulate the release of stored calcium ions within cells. This calcium ion release is an important regulatory process in cell signaling, apoptosis, cell proliferation, and other cellular physiological functions. Therefore, ITPRIP like 1 protein plays a significant role in the regulation of calcium ions within cells.

Research suggests that ITPRIP like 1 protein may be closely related to the normal functioning of the nervous system, muscle system, and other tissue organs. Further research will help better understand the mechanism of action of ITPRIP like 1 protein in biological processes, providing new insights and strategies for the treatment and prevention of related diseases.

Version:3.3 Revision

SDS-PAGE




On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

Bioactivity-ELISA

Bioactivity-ELISA

0.1 μg Biotinylated Cynomolgus ITPRIPL1 Protein; His-Avi Tag of per well



Biotinylated Cynomolgus ITPRIPL1 Protein; His-Avi Tag (Catalog # GM-87669RP) was immobilized at 1 $\mu\text{g}/\text{ml}$ (100 $\mu\text{L}/\text{well}$) on streptavidin precoated. Increasing concentrations of Anti-H_ITPRIPL1 mIgG1 Antibody(13B7A6H3) (Catalog # GM-87827AB) were added.