

Anti-EGFR hlgG1 Reference Antibody (Necbio)

Product Information

Product Name	Anti-EGFR hIgG1 Reference Antibody (Necbio)
Storage temp.	Store at 2-8°C short term (1-2 weeks). Store at \leq -20°C long term. Avoid repeated freeze-thaw.
Catalog# / Size	GM-87808MAB-1mg / 1 mg GM-87808MAB-5mg / 5 mg
	GM-87808MAB-25mg / 25 mg GM-87808MAB-50mg / 50 mg
	GM-87808MAB-100mg / 100 mg

Antibody Information

Expression System	СНО
Aggregation	< 5% as determined by SEC-HPLC
Purity	> 95% as determined by SDS-PAGE
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay
Sterility	0.2 µm Filtered
Target	EGFR
Clone	Necitumumab
Alternative Names	ERBB, ERBB1, ERRP, HER1, NISBD2, PIG61, mENA
Source/Isotype	Human IgG1(REEM),Kappa
Application	1

The epidermal growth factor receptor (EGFR; ErbB-1; HER1 in humans) is a transmembrane protein that is a receptor for members of the epidermal growth factor family (EGF family) of extracellular protein ligands. The epidermal growth factor receptor is a member of the ErbB family of receptors, a subfamily of four closely related receptor tyrosine kinases: EGFR (ErbB-1), HER2/neu (ErbB-2), Her 3 (ErbB-3) and Her 4 (ErbB-4). In many cancer types, mutations affecting EGFR expression or activity could result in cancer. Deficient signaling of the EGFR and other receptor tyrosine kinases in humans is associated with diseases such as Alzheimer's, while over-expression is associated with the development of a wide variety of tumors. Interruption of EGFR signalling, either by blocking EGFR binding sites on the extracellular domain of the receptor or by inhibiting intracellular tyrosine kinase activity, can prevent the growth of EGFR-expressing tumours and improve the patient's condition.

Formulation

Description


phosphate-buffered solution, pH 7.4.

Version:3.1




Data Examples

SDS-PAGE



On SDS-PAGE under reducing (R)/non-reducing(N-R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

SEC-HPLC



The purity of this product is more than 95% verified by SEC-HPLC

Version:3.1