

Mouse anti Human IgG1-DXD(Dar8)

Product Information


Product Name	Mouse anti Human IgG1-DXD(Dar8)
Storage temp.	Store at +4°C short term. Store at -80°C for 36 months, avoiding freeze/thaw cycles.
Catalog# / Size	GM-88039AB-100 / 100 µg GM-88039AB-1000 / 1 mg

Antibody Information

Expression System	CHO
Aggregation	< 5% as determined by SEC-HPLC
Purity	> 95% as determined by SDS-PAGE
Drug-to-Antibody Ratio (DAR)	8.0
Endotoxin	<1 EU/mg, determined by LAL gel clotting assay
Sterility	0.22 µm Filtered
Target	/
Clone	/
Alternative Names	/
Source/Isotype	Monoclonal Mouse IgG1 ,kappa
Application	ELISA, FACS, Cell-based Assay, Animal Model
Description	/
Formulation	Phosphate-buffered solution, pH 7.2-7.4.


Data Examples

SDS PAGE




After SDS-PAGE under reducing conditions, the gel was stained with Coomassie Blue overnight. The purity of the protein is greater than 95%.

SEC-HPLC



The purity of this product is greater than 95% verified by SEC-HPLC.

RP-HPLC



The drug-to-antibody ratio (DAR) of this product is 7.59.