

## **Product Sheet**

## Cynomolgus\_GPRC5D CHO-K1 Cell Line

Catalog number: GM-C40109

Version 3.3.1.250630

	Cynomolgus_GPRC5D CHO-K1 Cell Line is a clonal stable CHO-K1 cell line that	
Description	constitutively expresses the cynomolgus GPRC5D gene, constructed using lentiviral	
	technology.	
Quantity	5E6 Cells per vial,1 mL	
Product Format	1 vial of frozen cells	
Shipping	Shipped on dry ice	
Storage Conditions	Liquid nitrogen immediately upon receipt	
Target	Cynomolgus_GPRC5D	
Gene ID/Uniprot ID	XP_065379847.1	
Host Cell	CHO-K1	
<b>Recovery Medium</b>	F12K+10% FBS+1% P.S	
Growth medium	F12K+10% FBS+1% P.S+4 µg/mL Puromycin	
Note	None	
Freezing Medium	90% FBS+10% DMSO	
Growth properties	Adherent	
Growth Conditions	37°C, 5% CO <sub>2</sub>	
Growth Conditions	37°C, 5% CO <sub>2</sub>	
Growth Conditions Mycoplasma Testing	37°C, 5% CO <sub>2</sub> The cell line has been screened to confirm the absence of Mycoplasma species.	
Growth Conditions Mycoplasma Testing Safety considerations	37°C, 5% CO <sub>2</sub> The cell line has been screened to confirm the absence of Mycoplasma species. Biosafety Level 2	
Growth Conditions Mycoplasma Testing Safety considerations Note	37°C, 5% CO <sub>2</sub> The cell line has been screened to confirm the absence of Mycoplasma species. Biosafety Level 2 It is recommended to expand the cell culture and store a minimum of 10 vials at an early	



Genomeditech (Shanghai) Co.,Ltd. Order: +86 021-68455258/50432826/50432825 Toll-free: +86 400 627 9288 Email: service@genomeditech.com

## **Materials**

Reagent	Manufacturer/Catalogue No.
F12K	BOSTER/PYG0036
Fetal Bovine Serum	ExCell/FSP500
Pen/Strep	Thermo/15140-122
Puromycin	Genomeditech/GM-040401
Anti-GPRC5D Antibody (APC), Mouse Monoclonal	SinoBiological/24447-MM19-A

## Figures



Figure 1 | Cynomolgus\_GPRC5D CHO-K1 Cell Line (Cat. GM-C40109) was determined by flow cytometry using Anti-GPRC5D Antibody (APC), Mouse Monoclonal (SinoBiological/24447-MM19-A).

## **Cell Recovery**

#### Recovery Medium: F12K+10% FBS+1% P.S

To insure the highest level of viability, thaw the vial and initiate the culture as soon as possible upon receipt. If upon arrival, continued storage of the frozen culture is necessary, it should be stored in liquid nitrogen vapor phase and not at  $-70^{\circ}$ C. Storage at  $-70^{\circ}$ C will result in loss of viability.

- a) Thaw the vial by gentle agitation in a 37°C water bath. To reduce the possibility of contamination, keep the O-ring and cap out of the water. Thawing should be rapid (approximately 2 3 minutes).
- b) Remove the vial from the water bath as soon as the contents are thawed, and decontaminate by dipping in or spraying with 70% ethanol. All of the operations from this point on should be carried out under strict aseptic conditions.
- c) Transfer the vial contents to a centrifuge tube containing 5.0 mL complete culture medium and spin at approximately 176 x g for 5 minutes. Discard supernatant.
- d) Resuspend cell pellet with the recommended recovery medium. And dispense into appropriate culture dishes.

上海市浦东新区康威路 299 号 1 幢东区 505-507 邮编 201315 505-507,5th Floor, East District, Building 1,No.299 Kangwei Road, Pudong New Area, Shanghai 本公司产品仅供科研用途,严禁用于人体治疗! For research use only!



e) Incubate the culture at 37°C in a suitable incubator. A 5% CO<sub>2</sub> in air atmosphere is recommended if using the medium described on this product sheet.

## **Cell Freezing**

Freezing Medium: 90% FBS+10% DMSO

- a) Centrifuge at 176 x g for 3 minutes to collect cells.
- b) Resuspend the cells in pre-cooled freezing medium and adjust the cell density to 5E6 cells/mL.
- c) Aliquot 1 mL into each vial.
- d) Place the vial in a controlled-rate freezing container and store at -80°C for at least 1 day, then transfer to liquid nitrogen as soon as possible.

## Cell passage

Growth medium: F12K+10% FBS+1% P.S+4 µg/mL Puromycin

For the first 1 to 2 passages post-resuscitation, use the recovery medium. Once the cells have stabilized, switch to a growth medium.

- a) Remove and discard culture medium.
- b) Briefly rinse the cell layer with PBS to remove all traces of serum that contains trypsin inhibitor.
- c) Add 1.0 mL of 0.25% (w/v) Trypsin-EDTA solution to dish and observe cells under an inverted microscope until cell layer is dispersed (usually within 2 to 3 minutes at 37°C).
- Note: To avoid clumping do not agitate the cells by hitting or shaking the flask while waiting for the cells to detach.
   Cells that are difficult to detach may be placed at 37°C to facilitate dispersal.
- e) Add 2.0 mL of growth medium to mix well and aspirate cells by gently pipetting.
- f) After centrifugation, resuspend the pellet and add appropriate aliquots of the cell suspension to new culture vessels.
- g) Incubate cultures at 37°C.

#### Subcultivation Ratio: A subcultivation ratio of 1:4 - 1:5 is recommended

Medium Renewal: Every 2 to 3 days

## Notes

a) After the stabilization of the cell condition, there will be fewer dead cells post-passage, the cell growth rate will tend to stabilize, cell morphology will become uniform, and the cells will appear robust.

## Sequence

#### GPRC5D XP\_065379847.1

MYKDCIESTGDYFLPCDSEGPWGIVLESLAILGIVVTILLLLAFLFLMRKIQDCSQWNVLPTQLLFLLSVLGLF GLAFAFIIQLNQQTAPVRYFLFGVLFALCFSCLLAHASNLVKLVRGRVSFSWTTILCIAIGCSLLQVIIAIEYVTL IMTRGMMFVHMTPYQLNVDFVVLLVYVLFLMALTFFVSKATFCGPCENWKQHGRLIFITVLFSIIIWVVWIS  $\label{eq:mllrgnpqpqwddpvvcialvtnawvflllyivpelcilyrscrqecpsqghacpvtayqrsfqvenqelsrardsdgaeedvaltsfgtpiqpqtvdptqecfipraklspqqdagv*$ 

## **Related Products**

GPRC5D		
H_GPRC5D CHO-K1 Cell Line	H_GPRC5D HEK-293 Cell Line	
H_GPRC5D MC38 Cell Line		
Anti-H_GPRC5D hIgG4 Antibody(Talquetamab)		

## **License Agreement:**

# By purchasing and using this cell line product, the user voluntarily agrees to accept and abide by the following policies:

- This cell line product is restricted to research use only and shall not be used for any commercial purposes.
- This product is strictly prohibited from being used in the diagnosis or treatment of human or animal diseases, and shall not be directly used in experiments involving humans.
- Users are not permitted to modify the cell line in any way, nor to share, distribute, sell, sublicense, or otherwise transfer the licensed materials or their progeny to other laboratories, departments, research institutes, hospitals, universities, biotechnology companies, or any other third parties, except for research activities outsourced on behalf of the licensee. If the product is intended to be transferred to a third party, used for commercial development, preclinical or clinical drug functional validation, commercial production testing, or any other applications beyond the scope of this statement, prior written permission must be obtained from Genomeditech (Shanghai) Co.,Ltd. For details, please contact Genomeditech (Shanghai) Co.,Ltd.